

Urodynamics 2023



Anything is Possible
with the Right Approach



Index

Urodynamic Equipment	
Hermes	1.1-1.4
SmartDyn	1.5-1.9
SmartDyn Compact	1.10-1.14
BioSmart	1.15-1.16
Flowmeters	
SmartFlow	2.1-2.2
SmartScale	2.3-2.4
Aros	2.5-2.6
OruFlow	2.7
Bladder Scan	
Pinit	3.1-3.2
Urodynamic Chair	
Examination Chairs	4.1-4.10
Consumables	
Cystometry Catheters	5.1-5.6
Rectal Balloon Catheters	5.8-5.10
Catheter Pack (cysto + rectal)	5.12
Urodynamic Accesories	5.14-5.15
Pump Tubing	5.17
Pressure Transducers	5.19-5.20
Endo-Urology	5.22-5.25
WiScope (cystoscope and Ureteroscope)	
WiScope Imaging System - cystoscope and ureteroscope	6.1-6.2
Thermochemotherapy Equipment	
Thermochemotherapy System	7.1-7.2
Famale Incontinence	
DryNet	8.1
Probes	8.2

Urodynamic Equipment

Urodynamic Equipment

HERMES

The fully integrated Urodynamic System

The integrated design of the Hermes is a highly efficient urodynamic unit for a low cost investment, which facilitates a fast learning curve. Creo Medical complies with the highest quality standards, the quality of our catheters and after sales service ensure accurate results.



Workstation

- ▶ Wireless Connection of the Flowscale and the EMG to Patient Unit.
- ▶ Availability of studies and Nomograms from the ICS.
- ▶ Hygienic - Easy to Clean, with waterproof keyboard, and special designed protection to avoid damage of flowscale. Autoclavable penset for UPP studies.
- ▶ Optional Anorectal Manometry and Videourodynamics.
- ▶ Optional Integration with the Hospital Information System.

Integration

- ▶ Easy integration with current Hospital lay-out.
- ▶ Possibility to be used by different departments (Urology, Gynaecology) to reduce the effective cost per patient.

Quality

- ▶ Creo Medical guarantees the Highest Quality Standards in all its products.
- ▶ ISO 13845, CE Mark.

code	description
98600000	Hermes

Hermes Includes:




Filling Pump

- User friendly software control.
- Range from 0ml/min to 120ml/min.
- Automatic priming of the pump tubing set.




Filling Sensor

- Reports exact volume infused:
- Compatible with:
 - Different contrast media.
 - Catheters with different lumen sizes.
 - Avoids pump infusion errors.



Database / Analysis

- SQL Database engine ensures fast access to all studies:
- Complete analysis as defined by the ICS.
- User designed protocols and reports.
- NDL Technology (No Data Loss).



Nomograms / Reports


- Full range of Nomograms for adults and paediatric studies.
- Customised investigations and reports.



Network and HIS

- Wireless or cable connection to Hospital Network.
- DICOM compatible - Storage & Worklist.

Kit, Urodynamic pressure: (mandatory)



00002806

Pressure Transducers

- 3 colour coded pressure transducers.
- Urodynamic set for convenient flushing of the catheters.

FlowScale (mandatory)

Flow sensor - wireless (optional) or wired, specifically designed & sensitive for low flow-rates (ie paediatric or obstructed flow).



00004442

FlowScale with cable



98800000

FlowScale wireless

Optional Accessories for Hermes:



00004326

Puller

- Flexible, robust and easy to adapt to any position.
- Removable and sterilisable catheter holder to prevent cross contamination.

Adjustable Pole

- Secure and handy fixing system for catheters and pump tubing.
 - Quick and easy adjustment.
 - Hygienic.



00004616

EMG

- High Quality Signal.
- Wireless system allows free patient movement.
- Low electronic noise.
- Rechargeable.



00004408

00002404

Micturition Chair

- Height adjustable, folding arms, detachable seat and cover, removable funnel.


Micturition stand

- For male voiding studies, height adjustable, removable funnel.



98700000

Remote control




98500215

98500215D

Video-urodynamics

- High Resolution Processing.
- Available for Analogic (98500215) or Digital (98500215D) signals.
- Real time recording and synchronization.



98000003

Anorectal Manometry

- 4 or 5 Pressure channels.Vector Volume.
- Integrated air pressure pump for accurate perfusion.
- Foot pedal control.
- Configurable studies (RAIR, Pressure Profile,
- Fatigue Index, Sensitivity).

SmartDyn

Wireless Urodynamics System

The integrated design of the **SmartDyn** allows a fast learning curve. It is a high efficient device for a low cost investment.

Creo Medical after sales service and the quality of our catheters ensure the results of the measurements.

Creo Medical complies with the highest quality market standards.



Workstation

- ▶ Wireless Connection of the Flowscale and the EMG to Patient Unit.
- ▶ Availability of studies and Nomograms from the ICS.
- ▶ Hygienic - Easy to Clean, with waterproof keyboard, and special designed protection to avoid damage of flowscale. Autoclavable penset for UPP studies.
- ▶ Optional Anorectal Manometry and Videourodynamics.
- ▶ Optional Integration with the Hospital Information System.

Integration

- ▶ Easy integration with current Hospital lay-out.
- ▶ Possibility to be used by different departments (Urology, Gynaecology) to reduce the effective cost per patient.

Quality

- ▶ Creo Medical guarantees the Highest Quality Standards in all its products.
- ▶ ISO 13845, CE Mark.

code	description
98500000	SmartDyn

SmartDyn Includes:




Filling Pump

- User friendly software control.
- Range from 0ml/min to 120ml/min.
- Automatic priming of the pump tubing set.




Filling Sensor

- Reports exact volume infused:
- Compatible with:
 - Different contrast media.
 - Catheters with different lumen sizes.
 - Avoids pump infusion errors.



Database / Analysis

- SQL Database engine ensures fast access to all studies:
- Complete analysis as defined by the ICS.
- User designed protocols and reports.
- NDL Technology (No Data Loss).



Nomograms / Reports


- Full range of Nomograms for adults and paediatric studies.
- Customised investigations and reports.



Network and HIS

- Wireless or cable connection to Hospital Network.
- DICOM compatible - Storage & Worklist.

Kit, Urodynamic pressure: (mandatory)



00002806

Pressure Transducers

- 3 colour coded pressure transducers.
- Urodynamic set for convenient flushing of the catheters.

FlowScale (mandatory)

Flow sensor - wireless (optional) or wired, specifically designed & sensitive for low flow-rates (ie paediatric or obstructed flow).



00004442

FlowScale with cable



98800000

FlowScale wireless

Optional Accessories for SmartDyn:



00004326

Puller

- Flexible, robust and easy to adapt to any position.
- Removable and sterilisable catheter holder to prevent cross contamination.

Adjustable Pole

- Secure and handy fixing system for catheters and pump tubing.
 - Quick and easy adjustment.
 - Hygienic.



00004616

EMG

- High Quality Signal.
- Wireless system allows free patient movement.
- Low electronic noise.
- Rechargeable.



00004408

00002404

Micturition Chair

- Height adjustable, folding arms, detachable seat and cover, removable funnel.


Micturition stand

- For male voiding studies, height adjustable, removable funnel.



98700000

Remote control




98500215

98500215D

Video-urodynamics

- High Resolution Processing.
- Available for Analogic (98500215) or Digital (98500215D) signals.
- Real time recording and synchronization.



98000003

Anorectal Manometry

- 4 or 5 Pressure channels.Vector Volume.
- Integrated air pressure pump for accurate perfusion.
- Foot pedal control.
- Configurable studies (RAIR, Pressure Profile,
- Fatigue Index, Sensitivity).

Optional Accessories for SmartDyn:



98500214


Trolley workstation

- Can be used with SmartDyn to mount desktop or laptop PC, monitor, printer, keyboard etc. Matching colour with SmartDyn



00004393

Printer



00004325

All in one

SmartDyn Compact

Wireless Urodynamics System

The integrated design of the **SmartDyn Compact** allows a fast learning curve. It is a high efficient device for a low cost investment.

Creo Medical after sales service and the quality of our catheters ensure the results of the measurements.

Creo Medical complies with the highest quality market standards.



Workstation

- ▶ Wireless Connection of the Flowscale and the EMG to Patient Unit.
- ▶ Availability of studies and Nomograms from the ICS.
- ▶ Hygienic - Easy to Clean, with waterproof keyboard, and special designed protection to avoid damage of flowscale. Autoclavable penset for UPP studies.
- ▶ Optional Anorectal Manometry and Videourodynamics.
- ▶ Optional Integration with the Hospital Information System.

Integration

- ▶ Easy integration with current Hospital lay-out.
- ▶ Possibility to be used by different departments (Urology, Gynaecology) to reduce the effective cost per patient.

Quality

- ▶ Creo Medical guarantees the Highest Quality Standards in all its products.
- ▶ ISO 13845, CE Mark.

code	description
985C0000	SmartDyn Compact

SmartDyn Compact Includes:




Filling Pump

- User friendly software control.
- Range from 0ml/min to 120ml/min.
- Automatic priming of the pump tubing set.




Filling Sensor

- Reports exact volume infused:
- Compatible with:
 - Different contrast media.
 - Catheters with different lumen sizes.
 - Avoids pump infusion errors.



Database / Analysis

- SQL Database engine ensures fast access to all studies:
- Complete analysis as defined by the ICS.
- User designed protocols and reports.
- NDL Technology (No Data Loss).



Nomograms / Reports


- Full range of Nomograms for adults and paediatric studies.
- Customised investigations and reports.



Network and HIS

- Wireless or cable connection to Hospital Network.
- DICOM compatible - Storage & Worklist.

Kit, Urodynamic pressure: (mandatory)



00002806

Pressure Transducers

- 3 colour coded pressure transducers.
- Urodynamic set for convenient flushing of the catheters.

FlowScale (mandatory)

Flow sensor - wireless (optional) or wired, specifically designed & sensitive for low flow-rates (ie paediatric or obstructed flow).



00004442

FlowScale with cable



98800000

FlowScale wireless

Optional Accessories for SmartDyn Compact:



00004326

Puller

- Flexible, robust and easy to adapt to any position.
- Removable and sterilisable catheter holder to prevent cross contamination.

Adjustable Pole

- Secure and handy fixing system for catheters and pump tubing.
- Quick and easy adjustment.
- Hygienic.



00004616

EMG

- High Quality Signal.
- Wireless system allows free patient movement.
- Low electronic noise.
- Rechargeable.



00004408

00002404

Micturition Chair

- Height adjustable, folding arms, detachable seat and cover, removable funnel.


Micturition stand

- For male voiding studies, height adjustable, removable funnel.



98700000

Remote control




98500215

98500215D

Video-urodynamics

- High Resolution Processing.
- Available for Analogic (98500215) or Digital (98500215D) signals.
- Real time recording and synchronization.



98000003

Anorectal Manometry

- 4 or 5 Pressure channels.Vector Volume.
- Integrated air pressure pump for accurate perfusion.
- Foot pedal control.
- Configurable studies (RAIR, Pressure Profile,
- Fatigue Index, Sensitivity).

Optional Accessories for SmartDyn Compact:



98500214


Trolley workstation

- Can be used with SmartDyn to mount desktop or laptop PC, monitor, printer, keyboard etc. Matching colour with SmartDyn



00004393

Printer



00004325

All in one

BioSmart


Wireless Biofeedback Unit

- Wireless Communication Improves Ergonomics
- Patient Management
- Questionnaires
- User Target Configurations
- Standard and Customized Protocols
- Continuous Improvement Analysis
- Pelvic Floor Rehabilitation - Periform
- Biofeedback for Anal Sphincter - Anuform
- User Configurable Reports

Technical Specifications

EMG Input Channel
Pressure Input Channel
(only with SmartDyn Urodynamic System)

- ▶ Patient input protection: ~6kOhms
- ▶ Input range: 1mV full scale
- ▶ Input resolution: >100nV (averaged)
- Acceptable DC voltage: 600mV maximum
- Input impedance: 1010 ohms // 1pF
- Input impedance imbalance: >40%
- Common mode rejection: >100dB
- ▶ Bandwidth: 10Hz~1000Hz
- ▶ Optional Notch: 50Hz
- Input referred noise level: 2μV (averaged)
- Sampling rate: 2000 samples/second
- Resolution: 24bits (bit weight <30nV)



code	Description
00004821	BioSmart - Biofeedback Equipment

BioSmart Includes:



EMG

- High Quality Signal.
- Wireless system allows free patient movement.
- Low electronic noise.
- Rechargeable.




Vaginal probe



Anal probe



Middle connectors



Software Biofeedback

Flowmeters

Flowmeters

SmartFlow

Wireless Urodynamics

- The Smartflow has been specially designed as user friendly device, detecting the start of the miction and printing results automatically.
- The electronics used in the Smartflow allows the recording of very small flows, making easier the analysis with prostatic and paediatric patients.
- Configurable Report - MS Word.
- Patients Database.



Privacy:

- ▶ SmartFlow is a wireless urine flowmeter which allows privacy for patients improving the quality of the study.
- ▶ 10 to 30 meters signal depending on enviroment.

Control:

- ▶ Real time monitoring through a touchpanel LCD screen. All parameters according to international standards - ICS.








PC Management:

- ▶ Fiber optic connection with the computer.
- ▶ Nomograms - Liverpool, Siroky, Paediatric.
- ▶ Software 100% compatible with SmartDyn and Phoenix Plus.
- ▶ Configurable Report - MS Word.
- ▶ Patients Database.



code	description
00004440	SmartFlow

Optional Accesories and Consumables for SmartFlow

 <p>00004408 00002404</p>	<p>Micturition Chair · Height adjustable, folding arms, detachable seat and cover, removable funnel.</p> <p>Micturition stand · For male voiding studies, height adjustable, removable funnel.</p>
 <p>00004410</p>	<p>SmartFlow transportation case</p>
 <p>00004411</p>	<p>SmartFlow paper</p>
 <p>00004412</p>	<p>2L Jug for Smartflow</p>
 <p>00004415</p>	<p>500gr Calibration weight</p>
 <p>00004400 00004585B</p>	<p>Funnel for micturition chair</p> <p>Funnel for micturition stand</p>
 <p>00004413</p>	<p>Plasticover SmartFlow</p>

SmartScale

Wireless Flowmeter

- The SmartScale has been specially designed as user friendly device, detecting the start of the miction and printing results automatically.
- The electronics used in the Smartscale allows the recording of very small flows, making easier the analysis with prostatic and paediatric patients.
- Configurable Report - MS Word.
- Patients Database.

Privacy:

- ▶ **Smartscale** is a wireless urine flowmeter which allows privacy for patients improving the quality of the study.







PC Management:

- ▶ Fiber optic connection with the computer.
- ▶ Nomograms - Liverpool, Siroky, Paediatric.
- ▶ Software 100% compatible with SmartDyn and Phoenix Plus.
- ▶ Configurable Report - MS Word.
- ▶ Patients Database.



code	description
98800002	SmartScale

Optional Accesories for SmartScale

 <p>00004408 00002404</p>	<p>Micturition Chair · Height adjustable, folding arms, detachable seat and cover, removable funnel.</p> <p>Micturition stand · For male voiding studies, height adjustable, removable funnel.</p>
 <p>00004412</p>	<p>2L Jug for SmartScale</p>
 <p>00004415</p>	<p>500gr Calibration weight</p>
 <p>00004413</p>	<p>Plasticover SmartFlow</p>
 <p>98600062</p>	<p>Xbee</p>
 <p>00004400 00004585B</p>	<p>Funnel for micturition chair</p> <p>Funnel for micturition stand</p>

Aros

Combined flowmeter and EMG biofeedback system, for pelvic floor assessment.

- Suitable for the following studies:
 - Flowmeter - wireless
 - Liverpool, Siroky, Pediatrics
- High sensitivity electromyography
 - Surface electrodes, bipolar needles and rectal or vaginal cones
 - Wireless
- User configurable reports based on Microsoft Word. Unlimited number of templates.

Connectivity:

The Aros allows to work entirely with the computer history of the Hospital, using DICOM technology (optional). Allows DICOM-Worklist and DICOM-Store work modes.

Security (N.D.L):

The Aros has the most advanced technology N.D.L. which allows you to retrieve studies even if they have not been specifically recorded. The real-time processing of the data and its storage provide the user with security when carrying out the studies.

Database:

The Ring uses SQL as the database engine, a highly contrasted and versatile engine. This allows us to make manual and automatic backups in any location. In large databases, the AROS system has an intelligent search engine that will save time when reviewing studies.



code	description
98400000	Aros

001.005

Optional Accessories for Aros:

 <p>00004408 00002404</p>	<p>Micturition Chair · Height adjustable, folding arms, detachable seat and cover, removable funnel.</p> <p>Micturition stand · For male voiding studies, height adjustable, removable funnel.</p>
 <p>00004400 00004585B</p>	<p>Funnel for micturition chair</p> <p>Funnel for micturition stand</p>
 <p>00004412</p>	<p>2L Jug for Smartflow</p>
 <p>00004415</p>	<p>500gr Calibration weight</p>

UroFlow - Urinal

The integrated Uroflowmeter for men

- . Integrated team connected to the hospital's healthcare system
- . Automatic cleaning technology
- . High-precision transducer
- . Automatic Start/Stop
- . Touch control panel
- . Barcode identification
- . Remote monitoring via Smartphone or Tablet
- . Integration with the hospital's computer system - HIS



Benefits:

UroFlow saves Healthcare Professionals time because:

- There is no need to enter patient data.
- The cleaning process is completely automated.
- Patient monitoring is not required.
- Reports are automatically sent to the hospital's HIS (Hospital Information System).
- High-quality measurement:
- State-of-the-art technology equipment.
- Natural and private environment for the patient.
- No flow losses.

code	description
XXXX	UroFlow

Bladder Scan

Bladder Scan

Pinit

Non-invasive measurement of bladder volume

The Pinit Bladder Scanner is a "Real-time" portable ultrasound device dedicated to quick, non-invasive bladder volume measurement. The device is easy to operate by trained healthcare professionals.


Features

- ▶ Automatic volume measurement of bladder contents.
- ▶ Optional semi automatic mode - multipoint/trace contoured volume measurement – unique and excellent method for "difficult cases" and small volume bladders.
- ▶ Diagnostic scanning mode – classical 2D-mode, for urology (kidneys, bladder) and general abdominal examinations.
- ▶ Probe frequency: 2.7 MHz (Bladder Volume), 3.2 MHz (2D General)
- ▶ Quick real-time B-mode preview
- ▶ Touch-screen LCD control 5.7"
- ▶ Simple and intuitive user interface
- ▶ Mains power or battery power operated (internal rechargeable battery)
- ▶ Light and portable



11001000

Pinit ultrasound bladder scanner

 11001200	Bladder scan probe S-255 B
 11001001	Trolley/stand
 11001400	Pinit Paper roll pk/5
 11001007	Suitcase Pinit
11001300	Printer for urine volume measurement

Urodynamic Chairs

Urodynamic Chairs

Examination Chairs for Urology

Examination chair with electromechanical height adjustment, adjustment of pelvis section and back section.

Basic configuration

- ▶ Electromechanical movements controlled by easy operated cable foot switch
- ▶ Aluminium base and column cover powder coated, white (RAL 9010)
- ▶ Seat part and back rest frame made of polished stainless steel
- ▶ Floor levelling screws for secure standing
- ▶ Upholsteries with ergonomic form, 50 mm thickness
- ▶ Head and back upholstery prepared for neck upholstery with magnetic fixation
- ▶ Plastic covered backsides of upholstery parts
- ▶ Paper roll holder for paper rolls (up to 500 mm) at back rest
- ▶ Paper unreelable over back rest and seat part, or just over seat part
- ▶ Mounting device for side rails (25×10 mm), hand grips as well as a fixation device for leg holder system acc. to Goepel
- ▶ Extendable, inclinable, detachable white plastic basin
- ▶ Maximum permissible overall load: 200 kg
- ▶ Floor space required: 860 x 550 mm
- ▶ Power supply: 110/230 V / 50-60Hz , incl. 3 m cable



AM0400520

Examination Chair – Electrical Back,
Pelvis, Height

Examination Chairs for Urology

Examination chair with electromechanical height adjustment, adjustment of pelvis and back section, hand switch and memory function with three programmable positions for four operators.

Basic configuration

- ▶ Electromechanical height adjustment from 550 up to 900 mm
- ▶ Electromechanical adjustment of pelvis section up to 20°
- ▶ Electromechanical adjustment of back section up to 50°
- ▶ Electromechanical movements controlled by easy operated cable foot switch
- ▶ Additional hand switch for all electromechanical movements
- ▶ Memory function with three programmable positions for four operators
- ▶ Aluminium base and column cover powder coated, white (RAL 9010)
- ▶ Seat part and back rest frame made of polished stainless steel
- ▶ Floor levelling screws for secure standing
- ▶ Upholsteries with ergonomic form, 50 mm thickness
- ▶ Head and back upholstery prepared for neck upholstery with magnetic fixation
- ▶ Plastic covered backsides of upholstery parts
- ▶ Paper roll holder for paper rolls (up to 500 mm) at back rest
- ▶ Paper unreelable over back rest and seat part, or just over seat part
- ▶ Mounting device for side rails (25×10 mm), hand grips as well as a fixation device for leg holder system acc. to Goepel
- ▶ Extendable, inclinable, detachable white plastic basin
- ▶ Maximum permissible overall load: 200 kg
- ▶ Floor space required: 860 x 550 mm
- ▶ Power supply: 110/230 V / 50-60Hz , incl. 3 m cable



AM0400550

Examination Chair - Electrical Back,
Pelvis, Height & Memory

Examination Chairs for Urology and Gynaecology

Examination chair with electromechanical height adjustment, adjustment of pelvis and back section, leg holder adjustment, hand switch and memory function with three programmable positions for four operators.

Basic configuration

- ▶ Electromechanical movements controlled by easy operated cable foot switch
- ▶ Additional hand switch for all electromechanical movements
- ▶ Memory function with three programmable positions for four operators
- ▶ Aluminium base and column cover powder coated, white (RAL 9010)
- ▶ Seat part and back rest frame made of polished stainless steel
- ▶ Floor levelling screws for secure standing
- ▶ Upholsteries with ergonomic form, 50 mm thickness
- ▶ Head and back upholstery prepared for neck upholstery with magnetic fixation
- ▶ Plastic covered backsides of upholstery parts
- ▶ Paper roll holder for paper rolls (up to 500 mm) at back rest
- ▶ Paper unreelable over back rest and seat part, or just over seat part
- ▶ Mounting device for side rails (25×10 mm), hand grips as well as a fixation device for leg holder system acc. to Goepel
- ▶ Extendable, inclinable, detachable white plastic basin
- ▶ Maximum permissible overall load: 200 kg
- ▶ Floor space required: 860 x 550 mm
- ▶ Power supply: 110/230 V / 50-60Hz , incl. 3 m cable



AM0400570

Examination Chair – Electrical Back,
Pelvis, Height, Legs & Memory

Leg holder systems:



AM040508

Leg holder clamps

- White, RAL 9010
- For attachment at pelvis section
- Pair



AM081705

Leg holder according to Goepel Comfort

- With grey foamed PUR knee bolsters
- Pair



AM081717

Leg holders according to Goepel

- Made of stainless steel
- With integrated hand grip function and artificial leather cover
- For attachment at the chair's tothing elements
- Height adjustable
- Swivelling
- Grey foamed PUR footrests
- Stepless adjustment of the knee bolsters by ball joint
- Suitable for electric adjustment
- Pair



AM081726-2

Footrests

- Made of stainless steel
- With integrated hand grip function and artificial leather cover
- For attachment at the chair's tothing elements
- Height adjustable
- Swivelling
- Grey foamed PUR footrests
- Suitable for electric adjustment
- Pair



AM081726-1

Footrests Comfort

- Made of stainless steel
- With integrated hand grip function and artificial leather cover
- For attachment at the chair's tothing elements
- Stepless, gas spring supported one-hand height adjustment
- Pre-adjustment of the footrests by tothing elements at the chair
- Swivelling
- Grey foamed PUR footrests
- Pair

Leg holder systems:



AM081700

Leg holders according to Goepel Classic

- With black foamed PUR knee bolsters
- Pair
- Additional charge for artificial leather cover of knee bolsters (only 81700) in upholstery colour

Accessories for Gynaecology:



AM081047_200

Instrument drop


- For fixation at column of chair
- Extendable by means of telescopic rail (Tip-on)
- With stainless steel bowl (12.5 litres) and stainless steel basket for specula
- Dimensions (Stainless steel bowl): L x W x H = 300 x 235 x 200 mm
- Dimensions (Stainless steel basket): L x W x H = 230 x 175 x 140 mm



AM081051

Stainless steel board


- With extendible arm
- Swivelling
- Adaptable to the left or to the right
- Dimensions (L x W): 413 x 398 mm



AM081072

Gynaecological treatment set

- With extendible arm
- Swivelling
- Adaptable to the left or to the right
- Incl. additional accessories: 1 stainless steel storing-peel for instruments, 3 glass cups with lids, 3 glass bottles with lids, 4 instrument holster, 1 object slide




AM082347

Colposcope holder

- For systems Olympus, Kaps, Zeiss and Leisegang
- Adaptable to the left or to the right

Accessories for Urology:



AM091767

Drainage pan

- Made of stainless steel
- With tube



AM091765

Urodynamic seat cushion

- With central drainage channel



AM091762

Holder for uroflow system

- Incl. rack for chair attachment



AM092000

Holder for urodynamics patient unit

- For Albyn Medical Smartdyn"
- For side rail attachment



AM091766

Attachable upholsteries

- For enlargement of seat part for uro dynamic
- Dimensions (l x w): 170 x 150 mm



AM091608

Bucket

- Made of stainless steel
- Volume: 12 litres
- Incl. lid with tube-inlet



AM091617

Stainless steel disk

- With outlet and tube
- ø: 320 mm

Accessories for Urology:



AM091716

Hanging up rack

- Beside the chair for collecting basin or uroflow system
- Made of stainless steel



AM091706

Collecting basin

- Made of stainless steel
- With rinsing device
- Inflow electromagnetically controlled via foot switch
- $\varnothing = 320$ mm



AM093000

Elevators for flowmeter + funnel fastening system

AM001000

Remote control support on the left side of the patient

AM061511

Clamp for fixation of accessories with r(18mm), rotatable, made of stainless steel

AM081710

Extra cost of plastic lining (pair) for leg holder

AM400066

Extension for table transformation

AM400076

Extension for horizontal patient position

General accessories:



AM082349

Holder for examination luminaire

- Holder adaptable on the right side (Art.-no. 82349_1)
- Holder adaptable on the left side (Art.-no. 82349_2)
- Examination luminaire holder for examination chair with colposcope holder, adaptable at the right or at the left side (Art.-no. 82349_3)



AM082348

Examination luminaire

- LED technology (40'000 lx / 0.5 m)
- Illuminated field 180 mm
- Precise colour rendering Ra >93
- Light colours 3300 K, 3800 K or 4400 K
- Illumination strength dimmable in five stages
- Operating unit on light head
- Stepless positioning through flexible arm
- Slight thermal radiation in light



AM040070

Wall holder

- Made of stainless steel
- For attachment of footrests, leg holders, leg plates and further accessories to fix at side rail (25 x 10 mm)
- Dimensions (L x W x H): 300 x 235 x 200 mm



AM081805

Infusion rod

- For side rail attachment
- Stainless steel
- Incl. rotatable clamp



AM040023

Neck upholstery

- Incl. magnetic fixation



AM040220

Arm rest

- Swivelling
- Incl. artificial leather covered upholstery in colour of choice
- Dimensions (upholstery): L x W = 600 x 165 mm

General accessories:



AM061211

Arm rest

- Rotatable and tiltable by means of ball joint
- Height adjustable
- Horizontally slidable
- Upholstery electrically conductive
- Incl. clamp and fixture strap
- Dimensions (upholstery): L x W x H = 400 x 115 x 40 mm



AM028530

Hand grips

- 3 White RAL 9010
- Pair
- Additional charge for artificial leather cover of handgrips in colour of upholstery (Art.-no. 28531)



AM081043

Storing bowl

- Made of stainless steel
- Dimensions: (L x W x H): 325 x 265 x 65 mm



AM085223

Pelvis section side rails

- For accessory attachment
- Pair
- L= 320 mm



AM085222

Back section side rails

- For attachment of accessories
- Pair
- L= 600 mm



AM040020

PVC-protective cover

- For pelvis part upholstery



AM040046

Hand switch

- For control of height adjustment, pelvis section adjustment and back section adjustment
- Only for item 400520, cannot be retrofitted

General accessories:

 <p>AM040898</p>	<p>Mobility</p> <ul style="list-style-type: none"> · Consisting of 4 double leading castors · $\varnothing = 75\text{mm}$ · Lockable · Electrically conductive
 <p>AM041000</p>	<p>Powder coating of column</p> <ul style="list-style-type: none"> · Additional charge for powder coated cover of the column according to RAL colour chart
 <p>AM040000</p>	<p>Powder coating of base</p> <ul style="list-style-type: none"> · Additional charge for powder coated base according to RAL colour chart
 <p>AM042000</p>	<p>Column cover</p> <ul style="list-style-type: none"> · Made of polished stainless steel
 <p>AM040018</p>	<p>Doctor's foot rest</p> <ul style="list-style-type: none"> · Made of stainless steel
 <p>AM040678</p>	<p>Potential equalisation socket</p> <ul style="list-style-type: none"> · Cannot be retrofitted!
 <p>AM040677</p>	<p>Socket at the base</p> <ul style="list-style-type: none"> · 2 pieces · Cannot be retrofitted!

Consumables

Cystometry Catheters

Cystometry Catheters


Cystometry Catheters




Single lumen

Bladder fill

Used with code 00288114 (blow) for cystometry and pressure-flow studies. Tip is bent to facilitate tip of catheter 00288114, and allow easy separation once inside bladder. Higher filling rates are possible (>100ml/min) than with double-lumen catheters.



	description	material	diameter	hole from tip at	dimensions	QTY
00288010	Bladder fill	pebax	3.3mm/10Fr	15 mm	a=400mm	10





	description	material	diameter	hole from tip at	dimensions	QTY
95001001	Bladder fill	PVC	3.3mm/10Fr	15 mm	a=385mm	10
95001003	Bladder fill	PVC	2mm/6Fr	15 mm	a=400mm	10

Single lumen

Bladder pressure for flow study

An ideal catheter for pressure-flow studies, where minimum obstruction of the urethra during voiding is required. The catheter is long enough to connect directly to the Vesical pressure transducer. The tip fits inside the filling hold of code 00288019 (above) allowing easy introduction (catheters separate when inside bladder by "pulling-back")



	description	material	diameter	hole from tip at	dimensions	QTY
00288114	Bladder pressure	Pebax	1.7mm/5Fr	5/15 mm	a=1800mm	10
95001002	Bladder pressure	PVC	1.5mm/4.5Fr	10/20 mm	a=2000mm	25

Cystometry Catheters




Single lumen

Catheter set, bladder pressure pigtail and bladder fill

Tip of catheters should be able to fit inside hole of catheter code 00288020, and be separated by “pulling-back” once both catheters are inside the bladder.

Guidewire

- Material:** Guidewire: Steel Ø 0.5 mm
- Connections:** Guidewire: Male luer lock connector - colour CLEAR, to fit Albyn Medical female connector
- Packaging:** Pouch containing two catheters packaged “straight”; and guidewire. Guidewire located inside pressure catheter but not tightened, pigtail should not be straight.



	description	material	diameter	hole from tip at	dimensions	QTY
00288020	Bladder fill	pebax	1.7mm/5Fr	15 mm	a=1800mm	10
	Bladder pressure pigtail	50% pebax, 50% rilsan	3.3mm/10Fr	10/20/30 mm	a=450mm	



Cystometry Catheters



Double lumen




Double lumen catheters are ideal for filling/voiding cystometry studies, easy to introduce and comfortable for the patient. The above models should be used in conjunction with "Code 00288600 Connecting tube". The 6F model is ideal for Paediatric cystometry. Filling rates of up to 100ml/min are possible with both models.



A - tube length 60mm, to female luer lock connector (bladder pressure)
B - tube length 50mm, to female luer lock connector (fill)

	description	material	diameter	hole from tip at	dimensions	QTY
00288226	Cystometry Paediatric	Pebax	2mm/6Fr	A= 15mm B= 5mm	a= 400mm	10
00288228	Cystometry	Pebax	2.7mm/8Fr	A= 15mm B= 5mm	a= 400mm	10

Code 00288328 has an extended- pressure line, which allows connection directly to the pressure transducer, excluding the need for an additional adaptor/extension line. Filling rates up to 100ml/min possible.



A - tube length 1600mm, to female luer lock connector (bladder pressure)
B - tube length 200mm, to female luer lock connector (fill)

	description	material	diameter	hole from tip at	dimensions	QTY
00288326	Cystometry + ext line	Pebax	2mm/6Fr	A= 15mm B= 5mm	a= 400mm	10
00288328	Cystometry + ext line	Pebax	2.7mm/8Fr	A= 15mm B= 5mm	a= 400mm	10

Cystometry Catheters




Double lumen



A - tube length 1300mm, to female luer lock connector (bladder pressure)
 B - tube length 150mm, to female luer lock connector (fill)

	description	material	diameter	hole from tip at	dimensions	QTY
95012006	Cystometry + ext line	PVC	2mm/6Fr	A= 15mm B= 5mm	a= 700mm	10
95012008	Cystometry + ext line	PVC	2.7mm/8Fr	A= 15mm B= 5mm	a= 700mm	10
95012012	Cystometry + ext line	PVC	3.3mm/10Fr	A= 15mm B= 5mm	a= 700mm	10



A - tube length 1300mm, to female luer lock connector (bladder pressure)
 B - tube length 150mm, to female luer lock connector (fill)


	description	material	diameter	hole from tip at	dimensions	QTY
95013008	Cystometry Tiemann Tip	PVC	2.7mm/8Fr	A= 15mm B= 25mm	a= 470mm	10

Cystometry + UPP Catheter



Double lumen


This catheter may be used for cystometry and UPP studies. Following cystometry the filling lumen (blue) is disconnected from the filling pump and connected to the urethral pressure transducer using the extension line and "Y" piece (both included in pack). An infusion line is connected to the other port on the "Y" piece. Filling rates of up to 100ml/min possible. For Stress UPP measurements we recommend using a microtip catheter.



A - tube length 200mm, to female luer lock connector (bladder pressure)
 B - tube length 1400mm, to female luer lock connector (fill)

Extension line: For connection of catheter to urethral pressure transducer, length 1200mm

	description	material	diameter	hole from tip at	dimensions	QTY
00288426	Cystometry & Upp 2 lumen + ext line	Pebax	2mm/6Fr	A= 65mm B= 5mm	a= 400mm	10
00288428	Cystometry & Upp 2 lumen + ext line	Pebax	2.7mm/8Fr	A= 65mm B= 5mm	a= 400mm	10



A - 4 radial holes 90 degrees, Ø 0.95 mm and distal open
 B - 2 radial holes 180 degrees, Ø 0.95 mm

A - tube Ø 10.5 x 1.7 mm, Creo Medical winged female luer lock connector (colour RED), length 100 mm
 B - tube Ø 10.5 x 1.7 mm, Creo Medical winged female luer lock connector (colour BLUE), length 100 mm


	description	material	diameter	hole from tip at	dimensions	QTY
33288448	Cystometry catheter co-axial + UPP catheters 2 lumen + ext line	Pebax	2.7 mm/8 Fr	A= 65mm B= 5mm	a= 400mm	10

Urethrocystometry Catheter




Triple lumen


A single catheter for bladder filling and measuring Vesical and Urethral pressure, for filling/voiding cystometry and UPP's. Vesical (blue) and urethral (yellow) lumens connect to respective pressure transducers using code "00288600 Connecting tube".
Filling rates up to 100ml/min possible.



- A - tube length 60mm, to female luer lock connector (bladder pressure)
- B - tube length 50mm, to female luer lock connector (fill)
- C - tube length 70mm, to female luer lock connector, (urethral pressure)



	description	material	diameter	hole from tip at	dimensions	QTY
00288539	Urethrocystometry	Pebax	3mm/9Fr	A= 15mm B= 5mm C= 75mm	a= 400mm	10



- A - tube length 70mm, to female luer lock connector (bladder pressure)
- B - tube length 70mm, to female luer lock connector (fill)
- C - tube length 70mm, to female luer lock connector, (urethral pressure)

	description	material	diameter	hole from tip at	dimensions	QTY
00288639	Urethrocystometry	Polyurethane	3mm/9Fr	A= 10mm B= 70mm C= 70mm	a= 400mm	10

Rectal Balloon Catheters


Rectal Balloon Catheters

Rectal Balloon Catheters



Single lumen

Rectall balloon catheters provide the idea method for measuring abdominal pressure. The balloon prevents faeces from blocking the lumen, and provides good retention. Code 00288610 is suitable for paediatric studies.



Holes: 3 inside balloon to aid water filling/air removal


Connections: Female luer lock connector, colour red (rectal pressure) 1 way stopcock included

	description	material	diameter	balloon dimensions	dimensions	QTY
00288610	Rectal balloon	Tube:Pebax Balloon: Latex	4mm/ 12Fr	55mm Ø11mm	a= 180mm	10
00288611	Rectal balloon	Tube:Pebax Balloon: Latex	4mm/ 12Fr	70mm Ø15mm	a= 300mm	10

Holes: 3 inside balloon (distal, mid and proximal)

Connections: Female luer lock connector, colour red (rectal pressure)

	description	material	diameter	balloon dimensions	dimensions	QTY
00288612	Rectal balloon	Tube:Pebax Balloon: Silicone	4mm/ 12Fr	30mm Ø9mm	a= 180mm	10
00288613	Rectal balloon	Tube:Pebax Balloon: Silicone	4mm/ 12Fr	50mm Ø12mm	a= 300mm	10




	description	material	diameter	balloon dimensions	dimensions	QTY
95012004	Rectal balloon	Tube:PVC Balloon: PVC	4.5mm/ 13.5Fr	32mm Ø13,5mm	a= 2000mm	25

Rectal Balloon Catheters



Single lumen


The balloon on this catheter has a cut - to allow air to be flushed out during priming. The catheter is long enough to connect directly to the abdominal pressure transducer. The balloon is bonded to the tube to prevent it being detached and left inside the patient when the catheter is withdrawn.



Holes: 1 inside balloon

Connections: Female luer lock connector, colour red (rectal pressure)

	description	material	diameter	balloon dimensions	dimensions	QTY
95000000	Rectal balloon (cut)	Tube: PVC Balloon: PVC	1.7mm/ 5Fr	20mm Ø7mm	a= 1850mm	10
95000006	Rectal balloon (cut)	Tube: PVC Balloon: Silicone	1.7mm/ 5Fr	20mm Ø7mm	a= 1850mm	10





Introducer: 65mm length, 4.1 mm OD x 1.5mm ID

	description	material	diameter	balloon dimensions	dimensions	QTY
95012024	Rectal balloon (cut) Catheter with Introducer	Tube: PVC Balloon: PVC	4.5mm/ 13.5Fr	32mm Ø13,5mm	a= 2000mm	25

Rectal Balloon Catheters




Double lumen



A: 4.5Fr and 1300mm
B: 6Fr and 150mm

	description	material	diameter	holes from tip at	balloon dimensions	dimensions	QTY
95012015	Rectal balloon	Tube: PVC Balloon: PVC	3.3mm/ 10Fr	5/15mm	32mm Ø13,5mm	a= 700mm	10




A: 70mm
B: 70mm

	description	material	diameter	holes from tip at	balloon dimensions	dimensions	QTY
4114AE	Rectal balloon	Tube: Polyurethane Balloon: Polyisoprene	3.3mm/ 10Fr	10/20mm	45mm Ø11mm	a= 400mm	10

Catheter Pack (cysto + rectal)


Catheter Pack (cysto + rectal)

Catheter Pack




description QTY

95002006	Catheter Pack: · Double Lumen Cystometry Catheter (95012006) · Rectal PVC Balloon Catheter (95012004)	10
----------	---	----




description QTY

95002008	Catheter Pack: · Double Lumen Cystometry Catheter (95012008) · Rectal PVC Balloon Catheter (95012004)	10
----------	---	----



description QTY

95012010	Catheter Pack: · Double Lumen Cystometry Catheter (95012008) · Rectal PVC Balloon Catheter with introducer (95012024)	10
----------	---	----



description QTY

95012020	Catheter Pack: · Double Lumen Cystometry Catheter (95012006) · Rectal PVC Balloon Catheter with introducer (95012024)	10
----------	---	----

Urodynamic Accesories


Urodynamic Accesories

Urodynamic Accesories




One Way Valve Connection

Valve: Anti-return valve Ref: D703
 Connectors: Female Luer lock connector - Colour Transparent. Before valve.
 Male Luer lock connector - Colour Transparent. After valve.



	description	material	dimensions	QTY
9380000	One way Safety valve	PVC	120mm	20

The One way valve with extension is designed specifically to fit to pump tubing to avoid contamination of the PUMP by allowing circulation only in the filling diirection.





	description	material	dimensions	QTY
9380001	Infusion line extension with one-way valve	PVC	120mm	20

Connecting Tube

For use with Cystometry catheters 00288226 and 00288228, Urethrocystometry catheter 00288539, and rectal balloon catheters 00288610 and 00288611.

Connections: Female luer lock connector, colour clear
 Male luer lock connector, colour clear




	description	material	diameter	dimensions	QTY
00288600	Connecting Tube	PVC	2mm/6Fr	a= 1600mm	10

Urodynamic Accesories



Urethral Perfusion Line

- Material: $\varnothing a$ = PVC clear, 4mm external diameter
 $\varnothing b = \varnothing d = \varnothing e = \varnothing f$ = PVC clear, 2.5mm external diameter $\varnothing = 0.25\text{mm}$
- Connections: Connection at end of section $\varnothing a$ = spike (vented)
 Connection at end of section $\varnothing b$ = Creo Medical winged male luer lock connector - transparent
 Connection at end of section $\varnothing f$ = Creo Medical winged male luer lock connector - transparent
 red clamp fitted to section $\varnothing c$
 A = white cap
 The dimensions of the capillary tube should be in such way that for an applied pressure of 300 mmHg the flow through this tube must be 2ml/min.



	description	diameter	QTY
33288900	Urethral perfusion line	a = 1500mm b = 500mm c = 160mm d = 180mm	5

Pump Tubing

Pump Tubing


Pump Tubing Sets




Creo Medical

Both tubes are designed specifically for use with Albyn Medical Urodynamic system. The vented spike allows bottles to be used without the need for an additional vent or needle, and a clamp is fitted to help prevent unwanted spillage during priming. To help the user, the Green indicator on the luer connector matches with the Green connector on all catheters with a filling lumen. Code 00288710 will fit all systems, and the length of code 00288705 is tailored to fit "Phoenix Plus"

- Material: Sections "a" and "c" - PVC, clear
 Section "b" - Silicone rubber, clear
- Connections: Connection at end of section "a" = Spike (vented)
 Connection at end of section "c" = Rotating male luer lock colour clear,
 with green indicator at connector
 Clamp fitted to section "c"



	description	length	dimensions	QTY
00288705	Pump tube with clamp, for Phoenix Plus	3700mm	a = 500mm b = 200mm c = 3000mm	10
00288710	Pump tube with clamp, fo Creo Urodynamic systems	4200mm	a = 1000mm b = 200mm c = 3000mm	10



	description	length	dimensions	QTY
00288706	Set pump tube for Creo Medical equipment with pillow	3700mm	a = 500mm b = 200mm c = 2850mm d = 100mm	25

Pump tubing compatible with other brands also available, please contact us.

Pressure Transducers

Pressure Transducers

Pressure Transducers




Reusable pressure transducers and disposable domes provide the most cost-effective solution for urodynamic clinics

High quality reusable pressure transducer ideal for Urodynamics:

- 100mmHg test button
- Gold-plated non-corroding pins
- Silicone seals to prevent fluid from entering transducer-cable connection (with causes premature failure in other systems)

Packaging:

- Non-sterile, box containing 3, with colour-coded adhesive labels for Pves (blue), Pabd (red) and Pura (yellow).
- Transducer sensor protected with dummy dome



Pressure Transducer Domes

Disposable dome for reusable pressure transducer code 00002785
3-way stopcock fitted at one end, 1-way stopcock on other end




	description	QTY
00002787	Dome for pressure transducer	50



	description	QTY
00002787C	Dome for pressure transducer, Male/Male	50



Pressure Transducer for Anorrectal Module



	description	QTY
94122000	Pressure transducer for Anorrectal Module	50

Endo-Urology

UroAccess

Ureteral Access Sheath

4 configurations where the sheath is declined:

- In 2 diameters: 10 - 12Fr and 12 - 14Fr.
- In 2 lengths: 35 cm and 45 cm.

Characteristics

- ▶ **Optimised shape of the dilator tip** for easy insertion.
- ▶ **Bevelled access sheath end** for atraumatic transition between dilator and sheath



- ▶ Hydrophilic coating provides **excellent glide** during insertion.
- ▶ **Reinforced** and **radiopaque** sheath.
- ▶ Sheath is **graduated** along its length.

- ▶ **Luer-Lock connection** for contrast medium injection.
- ▶ **Locking system** of the expander in the sheath for a reliable coupling.



Reference	Ø Sheath interior (Fr)	Ø Sheath exterior (Fr)	Length (cm)	QTY
UAS-101235-041	10	12	35	4
UAS-101245-042	10	12	45	4
UAS-121435-038	12	14	35	4
UAS-121445-026	12	14	45	4

UroCatch


Extraction Basket

4 references **where the extraction basket is declined:**

- In **2 lengths: 90cm** and **120cm**.
- In **3 diameters: 1.9Fr, 2.2Fr** and **3.0Fr**.
- **With** or **without tip**.
- **Straight** or **helical** wire basket.

Characteristics

- ▶ Ergonomic handle, easy to **remove** and **re-assemble**.
- ▶ **360° rotatable** basket.
- ▶ **Radiopaque** sheath and basket.
- ▶ **Nitinol basket:** shape **memory material, resistant to deformation**.
- ▶ Models **without Tip** to **facilitate the capture of stones**, especially in the calcareous funds.



Reference	Ø Catheter (Fr)	Length (cm)	Basket	Ø Basket (mm)	Tip	QTY
USB-30090-001	3,0	90	4 helicoidal wires	13 - 16	✓	1
USB-30090-009	3,0	90	4 straight wires	13 - 16	✓	1
USB-19120-042	1,9	120	4 straight wires	11 - 16	✗	1
USB-22120-041	2,2	120	4 straight wires	11 - 16	✗	1

UroRig

Y-Tubing Irrigation Set

Characteristics

- ▶ Allows **good irrigation** during endoscopic procedures.
- ▶ **The manual pump with non-return valve allows:**
 - **Better flow control,**
 - Facilitates additional flow for **cleaning the surgical area.**
- ▶ **Two tips** to connect to **two irrigation bags.**
- ▶ **Three clamps** to **adapt the flow of fluids.**
- ▶ **T-connector** for **purging air** from inside the tubing.
- ▶ **Unitary** and **sterile** packaging.



Reference	Description	Length (cm)	Hand pump	QTY
UIS-225-121	Y-Tubing Irrigation Set	225	✓	20


UroGuide






Hydrophilic Guide Wire

- 5 available models:
- 2 diameters: 0.035" and 0.038"
 - 2 configurations: straight and curved
 - 1 reference straight at one end and curved at the other

Characteristics

- ▶ **Flexible tapered tip** for atraumatic operation.
- ▶ **Hydrophilic coating** for low friction insertion and advancement.
- ▶ **Nitinol core:**
 - **Flexible material, resistant to deformation.**
 - **Good radiopacity** providing good fluoroscopic visibility for controlled navigation.
- ▶ **Guide wire protective case with:**
 - **Luer-Lock connection** for easy and complete wetting of the guide wire before use.
 - **Introducer for centring and facilitating the insertion** of the guide wire into the inlet valve of the operating channel.



Reference	Ø (inch / mm)	Length (cm)	Length of flexible end	Shape of the flexible end	QTY
UGW-035150-005	0.035 / 0.89	150	1 x 3 cm		5
UGW-038150-006	0.038 / 0.97	150	1 x 3 cm		5
UGW-035150-011	0.035 / 0.89	150	1 x 3 cm		5
UGW-038150-012	0.038 / 0.97	150	1 x 3 cm		5
UGW-035150-040	0.035 / 0.89	150	2 x 5 cm		5

WiScope
(Cystoscope & Ureteroscope)

WiScope (Cystoscope & Ureteroscope)

WiScope Image System - Cystoscope & Ureteroscope


A high-performance, immediately available and 100% functional piece of equipment for your ureteroscopic and cystoscopic procedures.

In addition to its immediate availability, the single use ureteroscope and cystoscope WiScope allows:

- ▶ To quickly sequence interventions
- ▶ To manage emergencies
- ▶ To eliminate the risk of cross contamination
- ▶ To guarantee expected performance in terms of use and image quality




code	description
OM.1966.903	Trolley WiScope



code	description
00060001	WiScope Image System

Ureteroscope

- Support in both directions, 275° up, 275° down
- Minimum radius of curvature 8 mm
- Distal end Ø 7.4 Fr (2.47 mm)
- Outer sleeve Ø 8.6 Fr (2.85 mm)
- Max. outer cover Ø 9.5 Fr (3.25 mm)
- Working channel Ø 3.6 Fr (1.10 mm)
- Active length 670 mm
- Weight of the handle 80 gr
- Total weight 170 gr



Reference	Right-handed / Left-handed	Deflection	
00060000	Right-handed	Standard	Deflection lever upwards – distal tip « up » Deflection lever downwards – distal tip « down »
00060002	Left-handed	Standard	
00060003	Right-handed	Reverse	Deflection lever upwards – distal tip « down » Deflection lever downwards – distal tip « up »
00060004	Left-handed	Reverse	

Cystoscope

- Support in both directions, 180° up, 180° down
- Minimum radius of curvature 10 mm
- Distal end Ø 9.0 Fr (3.0 mm)
- Outer cover Ø 13.8 Fr (4.6 mm)
- Max. outer cover Ø 15 Fr (5.0 mm)
- Working channel Ø 6.6 Fr (2.2 mm)
- Active length 380 mm
- Weight of the handle 80 gr
- Total weight 170 gr



Reference	Right-handed / Left-handed	Deflection	
00060014	Right-handed	Standard	Deflection lever upwards – distal tip « up » Deflection lever downwards – distal tip « down »
00060013	Left-handed	Standard	
00060012	Right-handed	Reverse	Deflection lever upwards – distal tip « down » Deflection lever downwards – distal tip « up »
00060011	Left-handed	Reverse	

Thermochemotherapy Equipment

Thermochemotherapy Equipment

Thermochemotherapy System

Elmedical delivers a system (console and catheters) for minimally invasive non-ablative high temperature treatments (hyperthermia) for urological diseases.

BWT (Bladder Wall Thermochemotherapy) for Bladder cancer:

- ▶ BWT treats non-muscle invasive bladder cancer (NMIBC).
- ▶ BWT provides bladder wall thermochemotherapy with uniform hyperthermia all over the bladder by high flow of a heated Mitomycin-C solution.
- ▶ BWT provides conductive-heat which is easy to control within the bladder.
- ▶ BWT heat is administered precisely and uniformly throughout the bladder avoiding "hot spots" or ineffective "cold spots".

Bladder Wall Thermochemotherapy (BWT) treats non-muscle invasive bladder cancer (NMIBC) by combining uniform bladder wall hyperthermia with chemotherapy instillation.

Combination of heat and drug has synergistic effect:

- ▶ Heat selectively damages the cancer tissue (neoplastic tissue)
- ▶ Heat increases tissue drug uptake by 3 mechanisms:
 - Heat increases blood flow (more available to the tissue)
 - Heat increases cell membrane permeability to the drug
 - Heat accelerates drug metabolism (every 1 centigrade over 37 centigrade increase at 50%)

The procedure offering best patient's tolerance:

- ▶ No pain; performed with topical anesthesia only.
- ▶ No post-procedure Catheter.
- ▶ No rectal, bladder or tissue temperature probes.
- ▶ Patient returns back to normal activities the same day.



BWT inherent features:

- ▶ Better controlled conductive heat leads to better drug stability.
- ▶ The high rate circulation in the bladder continuously brings "fresh" heated drug in contact with all bladder wall.
- ▶ Uniform and precise bladder wall temperature maintained all over the bladder while there are:
 - No "hot spots" (burns in posterior bladder wall), and
 - No "cold spots" (ineffective low heating near the bladder dome or trigone).
 - The uniform and precise temperature enables increase of the targeted temperature to 44degC without adverse events.
 - The absence of burns leads to increased patient tolerability to the treatment.
 - Non-traumatic insertion of the 18Fr flexible catheter which has no rigid parts.

code	description
97100000	Thermochemotherapy System

Optional Accessories for Thermochemotherapy System




Trolley for Elmedical PelvixTT

97100003



97100001

Catheter set BWT 18Fr. Silicone catheter



97100002

DuoThermia catheter, for Prostitus treatment



97100004


UniThermia catheter set 16Fr. Tiemen

Female Incontinence

Famale Incontinence

DryNet System

DryNet is the ultimate solution for female stress urinary incontinence, a high quality, ultra-soft TOT mesh that minimises any post-surgical complications.



- ▶ Controlled elasticity
- ▶ Ergonomic cone-sleeve connection
- ▶ Recovers original shape after elongation
- ▶ Ultra soft to the touch

code	description
99100050	DryNet System

Periform Intravaginal Probe

The Periform Intravaginal Probe has been designed entirely for more comfortable and effective use of both clinical and home muscle stimulation systems. The improved interface also allows for better utilisation of EMG bio-reporting equipment and, thanks to its revolutionary pelvic floor contraction indicator, provides you and your medical consultant with valuable information to assist in the healthy development of the pelvic floor. Periform's revolutionary patented shape means that for most users the probe will not move out of place. This allows for more accurate assessment in a standing position and allows complete freedom of movement when using a muscle stimulator at home.



code

description

00900113

Periform Intravaginal Probe

Anuform Intraanal Probe

The Anuform Anal electrode can be used in conjunction with more EMS (electrical muscle stimulation) units and pelvic floor exercises to help with stress relief, urge, faecal and mixed incontinence for men and women. Also use it as a smaller alternative to vaginal insertion in conjunction with female incontinence relief units.



code

description

00900114

Anuform Intraanal Probe



Creo Medical authorised local distributor

Creo Medical S.L.

Pol. Ind. Cordovilla, Calle D. Nº1
31191 Cordovilla, Navarra (España)

Tel.: +34 948 293 324

E-mail: info.es@creomedical.com

Anything is Possible
with the Right Approach

CURO_E_3. Rev03. Date 11-04-23